



Human CADM4 blocking peptide (CDBP0658)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-CADM4 (C Terminus) antibody
Antigen Description	CADM4 (cell adhesion molecule 4) is a protein-coding gene. Diseases associated with CADM4 include renal clear cell carcinoma, and colon cancer. An important paralog of this gene is CADM3.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	CADM4 cell adhesion molecule 4 [Homo sapiens]
Official Symbol	CADM4
Synonyms	CADM4; cell adhesion molecule 4; IGSF4C, immunoglobulin superfamily, member 4C; Necl 4; nectin like 4; SynCAM4; TSLL2; TSLC1-like 2; nectin-like 4; TSLC1-like protein 2; nectin-like protein 4; immunoglobulin superfamily member 4C; immunoglobulin superfamily, member 4C; NECL4; IGSF4C; Necl-4; synCAM4;
Entrez Gene ID	199731

mRNA Refseq	NM_145296
Protein Refseq	NP_660339
UniProt ID	Q8NFZ8
Chromosome Location	19q13.32