



Human CBX1 blocking peptide (CDBP0708)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-CBX1/HP1-Beta antibody
Antigen Description	This gene encodes a highly conserved nonhistone protein, which is a member of the heterochromatin protein family. The protein is enriched in the heterochromatin and associated with centromeres. The protein has a single N-terminal chromodomain which can bind to histone proteins via methylated lysine residues, and a C-terminal chromo shadow-domain (CSD) which is responsible for the homodimerization and interaction with a number of chromatin-associated nonhistone proteins. The protein may play an important role in the epigenetic control of chromatin structure and gene expression. Several related pseudogenes are located on chromosomes 1, 3, and X. Multiple alternatively spliced variants, encoding the same protein, have been identified.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	CBX1 chromobox homolog 1 [Homo sapiens]
Official Symbol	CBX1

Synonyms	CBX1; chromobox homolog 1; chromobox homolog 1 (Drosophila HP1 beta); chromobox protein homolog 1; CBX; HP1 beta homolog (Drosophila); HP1 BETA; HP1Hs beta; M31; MOD1; HP1 beta homolog; modifier 1 protein; heterochromatin protein 1-beta; heterochromatin protein p25 beta; heterochromatin protein 1 homolog beta; chromobox homolog 1 (HP1 beta homolog Drosophila); p25beta; HP1-BETA; HP1Hsbeta; HP1Hs-beta;
Entrez Gene ID	10951
mRNA Refseq	NM_001127228
Protein Refseq	NP_001120700
UniProt ID	P83916
Chromosome Location	17q21.32
Function	chromatin binding; enzyme binding; identical protein binding; protein binding;