



Human CBX8 blocking peptide (CDBP0712)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-CBX8 antibody
Antigen Description	CBX8 (chromobox homolog 8) is a protein-coding gene. Diseases associated with CBX8 include childhood leukemia, and mouth disease. GO annotations related to this gene include ubiquitin-protein ligase activity and single-stranded RNA binding. An important paralog of this gene is CBX6.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	CBX8 chromobox homolog 8 [Homo sapiens]
Official Symbol	CBX8
Synonyms	CBX8; chromobox homolog 8; chromobox homolog 8 (Drosophila Pc class); chromobox protein homolog 8; HPC3; Pc class 3 homolog (Drosophila); PC3; polycomb 3; RC1; rectachrome 1; Pc class 3 homolog; chromobox homolog 8 (Pc class homolog, Drosophila);
Entrez Gene ID	57332

mRNA Refseq	NM_020649
Protein Refseq	NP_065700
UniProt ID	Q9HC52
Chromosome Location	17q25.3
Function	methylated histone residue binding; protein binding; single-stranded RNA binding; ubiquitin-protein ligase activity;