



Human CCDC106 blocking peptide (CDBP0715)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-CCDC106 antibody
Antigen Description	CCDC106 (coiled-coil domain containing 106) is a protein-coding gene.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	CCDC106 coiled-coil domain containing 106 [Homo sapiens]
Official Symbol	CCDC106
Synonyms	CCDC106; coiled-coil domain containing 106; coiled-coil domain-containing protein 106; HSU79303; protein predicted by clone 23882; ZNF581;
Entrez Gene ID	29903

mRNA Refseq	NM_013301
Protein Refseq	NP_037433
UniProt ID	Q9BWC9
Chromosome Location	19q13.42