



Human CDC37L1 blocking peptide (CDBP0740)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-CDC37L1/HARC antibody
Antigen Description	CDC37L1 is a cytoplasmic phosphoprotein that exists in complex with HSP90 (HSPCA; MIM 140571) as well as several other proteins involved in HSP90-mediated protein folding (Scholz et al., 2001 (PubMed 11413142)).
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	CDC37L1 cell division cycle 37 homolog (S. cerevisiae)-like 1 [Homo sapiens]
Official Symbol	CDC37L1
Synonyms	CDC37L1; cell division cycle 37 homolog (S. cerevisiae)-like 1; CDC37 cell division cycle 37 homolog (S. cerevisiae) like 1 , cell division cycle 37 homolog (S. cerevisiae) like 1; hsp90 co-chaperone Cdc37-like 1; CDC37B; FLJ20639; HARC; Hsp90-associating relative of Cdc37; CDC37 cell division cycle 37 homolog-like 1; RP11-6J24.5;
Entrez Gene ID	55664

mRNA Refseq	NM_017913
Protein Refseq	NP_060383
UniProt ID	Q7L3B6
Chromosome Location	9p24.1