



Human CHD7 blocking peptide (CDBP0784)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-CHD7 antibody
Antigen Description	This gene encodes a protein that contains several helicase family domains. Mutations in this gene have been found in some patients with the CHARGE syndrome.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	CHD7 chromodomain helicase DNA binding protein 7 [Homo sapiens]
Official Symbol	CHD7
Synonyms	CHD7; chromodomain helicase DNA binding protein 7; FLJ20357; FLJ20361; KIAA1416;
Entrez Gene ID	338420
UniProt ID	Q9P2D1
Chromosome Location	8q12.2