



Human CLEC16A blocking peptide (CDBP0818)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-CLEC16A (aa1040-53) antibody
Antigen Description	This gene encodes a member of the C-type lectin domain containing family. Single nucleotide polymorphisms in introns of this gene have been associated with diabetes mellitus, multiple sclerosis and rheumatoid arthritis. Multiple transcript variants encoding different isoforms have been found for this gene.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	CLEC16A C-type lectin domain family 16, member A [Homo sapiens]
Official Symbol	CLEC16A
Synonyms	CLEC16A; C-type lectin domain family 16, member A; KIAA0350; protein CLEC16A; Gop 1; Gop-1; MGC111457;
Entrez Gene ID	23274

mRNA Refseq	NM_001243403
Protein Refseq	NP_001230332
UniProt ID	Q2KHT3
Chromosome Location	16p13.13