



Human DDA1 blocking peptide (CDBP0982)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-DDA1 (aa32-42) antibody
Antigen Description	DDA1 (DET1 and DDB1 associated 1) is a protein-coding gene.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	DDA1 DET1 and DDB1 associated 1 [Homo sapiens]
Official Symbol	DDA1
Synonyms	DDA1; DET1 and DDB1 associated 1; C19orf58, chromosome 19 open reading frame 58; DET1- and DDB1-associated protein 1; MGC2594; PCIA1; PCIA-1; cross-immune reaction antigen PCIA1; placenta cross-immune reaction antigen 1; C19orf58;
Entrez Gene ID	79016
mRNA Refseq	NM_024050

Protein Refseq	NP_076955
UniProt ID	Q9BW61
Chromosome Location	19p13.11