



Human DEPDC4 blocking peptide (CDBP0994)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-DEPDC4 antibody
Antigen Description	DEPDC4 (DEP domain containing 4) is a protein-coding gene. An important paralog of this gene is DEPDC1.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	DEPDC4 DEP domain containing 4 [Homo sapiens]
Official Symbol	DEPDC4
Synonyms	DEPDC4; DEP domain containing 4; DEP domain-containing protein 4; DEP.4; FLJ33505;
Entrez Gene ID	120863
mRNA Refseq	NM_152317
Protein Refseq	NP_689530

UniProt ID	Q8N2C3
-------------------	--------

Chromosome Location	12q23
----------------------------	-------