



Human FRMPD3 blocking peptide (CDBP1298)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-FRMPD3 antibody
Antigen Description	FRMPD3 (FERM and PDZ domain containing 3) is a protein-coding gene. An important paralog of this gene is FRMPD1.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	FRMPD3 FERM and PDZ domain containing 3 [Homo sapiens]
Official Symbol	FRMPD3
Synonyms	FRMPD3; FERM and PDZ domain containing 3; FERM and PDZ domain-containing protein 3; KIAA1817; RP5 1070B1.1; RP5-1070B1.1;

Entrez Gene ID	84443
UniProt ID	Q5JV73
Chromosome Location	Xq22