



Human GPRIN3 blocking peptide (CDBP1415)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-GPRIN3 antibody
Antigen Description	GPRIN3 (GPRIN family member 3) is a protein-coding gene. Diseases associated with GPRIN3 include gastroesophageal reflux disease. An important paralog of this gene is GPRIN1.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	GPRIN3 GPRIN family member 3 [Homo sapiens]
Official Symbol	GPRIN3
Synonyms	GPRIN3; GPRIN family member 3; G protein-regulated inducer of neurite outgrowth 3; FLJ42625; GRIN3;
Entrez Gene ID	285513
mRNA Refseq	NM_198281

Protein Refseq	NP_938022
UniProt ID	Q6ZVF9
Chromosome Location	4q22.1