



# Human HS1BP3 blocking peptide (CDBP1516)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Blocking/Immunizing peptide for anti-HS1-binding protein 3 antibody
<b>Antigen Description</b>	The protein encoded by this gene shares similarity with mouse Hs1bp3, an Hcls1/Hs1-interacting protein that may be involved in lymphocyte activation. [provided by RefSeq, Jul 2008]
<b>Species</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Apuri, BL, ELISA
<b>Format</b>	Lyophilized powder
<b>Size</b>	100 µg
<b>Preservative</b>	None
<b>Storage</b>	Shipped at ambient temperature, store at -20°C.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">HS1BP3 HCLS1 binding protein 3 [ Homo sapiens ]</a>
<b>Official Symbol</b>	HS1BP3
<b>Synonyms</b>	HS1BP3; HCLS1 binding protein 3; HCLS1-binding protein 3; HS1 BP3; FLJ14249; HSP1BP-3; HS1-binding protein 3; ETM2; HS1-BP3;
<b>Entrez Gene ID</b>	<a href="#">64342</a>
<b>mRNA Refseq</b>	<a href="#">NM_022460</a>

<b>Protein Refseq</b>	<a href="#">NP_071905</a>
<b>UniProt ID</b>	Q53T59
<b>Chromosome Location</b>	2p24.1
<b>Function</b>	phosphatidylinositol binding;