



Mouse Ier5 blocking peptide (CDBP1551)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Ier5 (mouse) antibody
Antigen Description	This gene encodes a protein that is similar to other immediate early response proteins. In the mouse, a similar gene may play an important role in mediating the cellular response to mitogenic signals. Studies in rats found the expression of a similar gene to be increased after waking and sleep deprivation. [provided by RefSeq, Jul 2008]
Species	Mouse
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	IER5 immediate early response 5 [Homo sapiens]
Official Symbol	IER5
Synonyms	IER5; immediate early response 5;
Entrez Gene ID	30015