



Human JMJD4 blocking peptide (CDBP1648)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-JMJD4 antibody
Antigen Description	JMJD4 (jumonji domain containing 4) is a protein-coding gene. An important paralog of this gene is JMJD6.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	JMJD4 jumonji domain containing 4 [Homo sapiens]
Official Symbol	JMJD4
Synonyms	JMJD4; jumonji domain containing 4; jmjC domain-containing protein 4; FLJ12517; jumonji domain-containing protein 4; MGC129896;

Entrez Gene ID	65094
mRNA Refseq	NM_023007
Protein Refseq	NP_075383
UniProt ID	Q9H9V9
Chromosome Location	1q42.13