



Human LXN blocking peptide (CDBP1734)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Latexin antibody
Antigen Description	This gene encodes the only known protein inhibitor of zinc-dependent metallopeptidases. [provided by RefSeq, Oct 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	LXN latexin [Homo sapiens]
Official Symbol	LXN
Synonyms	LXN; latexin; MUM; tissue carboxypeptidase inhibitor; endogenous carboxypeptidase inhibitor; ECI; TCI;
Entrez Gene ID	56925
mRNA Refseq	NM_020169

Protein Refseq	NP_064554
UniProt ID	Q9BS40
Chromosome Location	3q25.32
Function	enzyme inhibitor activity; metalloendopeptidase inhibitor activity;