



Human LAX1 blocking peptide (CDBP1736)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-LAX1 antibody
Antigen Description	LAX1 (lymphocyte transmembrane adaptor 1) is a protein-coding gene. Diseases associated with LAX1 include blepharochalasis, and breast giant fibroadenoma. GO annotations related to this gene include SH2 domain binding and protein kinase binding.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	LAX1 lymphocyte transmembrane adaptor 1 [Homo sapiens]
Official Symbol	LAX1
Synonyms	LAX1; lymphocyte transmembrane adaptor 1; lymphocyte transmembrane adaptor 1; FLJ20340; LAT like membrane associated protein; LAX; linker for activation of x cells; LAT-like membrane associated protein; membrane-associated adapter protein LAX;
Entrez Gene ID	54900

mRNA Refseq	NM_001136190
Protein Refseq	NP_001129662
UniProt ID	Q8I WV1
Chromosome Location	1q32.1
Pathway	T Cell Receptor Signaling Pathway, organism-specific biosystem;
Function	SH2 domain binding; protein binding; protein kinase binding;