



Human NDFIP1 blocking peptide (CDBP1983)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-NDFIP1 antibody
Antigen Description	The protein encoded by this gene belongs to a small group of evolutionarily conserved proteins with three transmembrane domains. It is a potential target for ubiquitination by the Nedd4 family of proteins. This protein is thought to be part of a family of integral Golgi membrane proteins. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	NDFIP1 Nedd4 family interacting protein 1 [Homo sapiens]
Official Symbol	NDFIP1
Synonyms	NDFIP1; Nedd4 family interacting protein 1; NEDD4 family-interacting protein 1; MGC10924; N4WBP5; Nedd4 WW domain-binding protein 5; putative MAPK-activating protein PM13;

breast cancer-associated protein SGA-1M; putative NFKB and MAPK-activating protein;
putative NF-kappa-B-activating protein 164;

Entrez Gene ID	80762
mRNA Refseq	NM_030571
Protein Refseq	NP_085048
UniProt ID	Q9BT67
Chromosome Location	5q31.3
Function	WW domain binding; protein binding; signal transducer activity;