



Human GBAS blocking peptide (CDBP2047)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-NIPSNAP2 antibody
Antigen Description	This gene encodes a member of the NipSnap family of proteins that may be involved in vesicular transport. The encoded protein is localized to mitochondria and plays a role in oxidative phosphorylation. A pseudogene of this gene is located on the long arm of chromosome 2. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name [GBAS glioblastoma amplified sequence \[Homo sapiens \]](#)

Official Symbol	GBAS
Synonyms	GBAS; glioblastoma amplified sequence; protein NipSnap homolog 2; NIPSNAP2; 4-nitrophenylphosphatase domain and non-neuronal SNAP25-like 2;
Entrez Gene ID	2631
mRNA Refseq	NM_001202469
Protein Refseq	NP_001189398
UniProt ID	O75323
Chromosome Location	7p12
Function	protein binding;