



Human NIPSNAP3B blocking peptide (CDBP2050)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-NIPSNAP3B antibody
Antigen Description	NIPSNAP3B belongs to a family of proteins with putative roles in vesicular trafficking (Buechler et al., 2004 [PubMed 15177564]).[supplied by OMIM, Mar 2008]
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	NIPSNAP3B nipsnap homolog 3B (C. elegans) [Homo sapiens]
Official Symbol	NIPSNAP3B
Synonyms	NIPSNAP3B; nipsnap homolog 3B (C. elegans); protein NipSnap homolog 3B; FLJ11275;

SNAP1; NIPSNAP3;

Entrez Gene ID [55335](#)

mRNA Refseq [NM_018376](#)

Protein Refseq [NP_060846](#)

UniProt ID Q9BS92

Chromosome Location 9q31.3
