



Human ODZ3 blocking peptide (CDBP2113)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-ODZ3/Teneurin-3 antibody
Antigen Description	TENM3 (teneurin transmembrane protein 3) is a protein-coding gene. Diseases associated with TENM3 include coloboma, and microphthalmia. GO annotations related to this gene include protein heterodimerization activity and protein homodimerization activity. An important paralog of this gene is TENM2.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ODZ3 odz, odd Oz/ten-m homolog 3 (Drosophila) [Homo sapiens]
Official Symbol	ODZ3
Synonyms	ODZ3; odz, odd Oz/ten-m homolog 3 (Drosophila); teneurin-3; KIAA1455; Ten M3; ten-3; tenascin-M3; ODZ3-like protein; odd Oz/Ten-m homolog 3; protein Odd Oz/ten-m homolog 3; TNM3; Ten-m3; FLJ10474; FLJ10886;
Entrez Gene ID	55714

mRNA Refseq	NM_001080477
Protein Refseq	NP_001073946
UniProt ID	Q9P273
Chromosome Location	4q35