



Human ORAI3 blocking peptide (CDBP2127)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-ORAI3 antibody
Antigen Description	ORAI3 (ORAI calcium release-activated calcium modulator 3) is a protein-coding gene. Diseases associated with ORAI3 include estrogen-receptor negative breast cancer, and breast cancer. GO annotations related to this gene include protein binding. An important paralog of this gene is ORAI2.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	ORAI3 ORAI calcium release-activated calcium modulator 3 [Homo sapiens]
Official Symbol	ORAI3
Synonyms	ORAI3; ORAI calcium release-activated calcium modulator 3; TMEM142C, transmembrane

protein 142C; protein orai-3; MGC13024; transmembrane protein 142C; TMEM142C;

Entrez Gene ID	93129
mRNA Refseq	NM_152288
Protein Refseq	NP_689501
UniProt ID	Q9BRQ5
Chromosome Location	16p11.2
Function	protein binding;