



Human PCSK9 blocking peptide (CDBP2225)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-PCSK9 antibody
Antigen Description	This gene encodes a member of the subtilisin-like proprotein convertase family, which includes proteases that process protein and peptide precursors trafficking through regulated or constitutive branches of the secretory pathway. The encoded protein undergoes an autocatalytic processing event with its prosegment in the ER and is constitutively secreted as an inactive protease into the extracellular matrix and trans-Golgi network. It is expressed in liver, intestine and kidney tissues and escorts specific receptors for lysosomal degradation. It plays a role in cholesterol and fatty acid metabolism. Mutations in this gene have been associated with autosomal dominant familial hypercholesterolemia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2014]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	PCSK9 proprotein convertase subtilisin/kexin type 9 [Homo sapiens]
Official Symbol	PCSK9

Synonyms	PCSK9; proprotein convertase subtilisin/kexin type 9; HCHOLA3, hypercholesterolemia, autosomal dominant 3; FH3; NARC 1; NARC1; NARC-1; HCHOLA3;
Entrez Gene ID	353175
UniProt ID	Q8NBP7
Chromosome Location	1p34.1-p32