



Human TBKBP1 blocking peptide (CDBP2403)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-Prosapip2 antibody
Antigen Description	TBKBP1 is an adaptor protein that binds to TBK1 (MIM 604834) and is part of the interaction network in the TNF (MIM 191160)/NFKB (see MIM 164011) pathway (Bouwmeester et al., 2004 [PubMed 14743216]).[supplied by OMIM, Mar 2008]
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	TBKBP1 TBK1 binding protein 1 [Homo sapiens]
Official Symbol	TBKBP1
Synonyms	TBKBP1; TBK1 binding protein 1; TANK-binding kinase 1-binding protein 1; KIAA0775; ProSAPiP2; TBK1-binding protein 1; SINTBAD;

Entrez Gene ID	9755
mRNA Refseq	NM_014726
Protein Refseq	NP_055541
UniProt ID	A7MCY6
Chromosome Location	17q21.32
Pathway	RIG-I-like receptor signaling pathway, organism-specific biosystem; RIG-I-like receptor signaling pathway, conserved biosystem;