



Human SKA3 blocking peptide (CDBP2690)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-SKA3 antibody
Antigen Description	This gene encodes a component of the spindle and kinetochore-associated protein complex that regulates microtubule attachment to the kinetochores during mitosis. The encoded protein localizes to the outer kinetechore and may be required for normal chromosome segregation and cell division. Alternative splicing results in multiple transcript variants. (provided by RefSeq, Oct 2009)
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	SKA3 spindle and kinetochore associated complex subunit 3 [Homo sapiens]
Official Symbol	SKA3
Synonyms	SKA3; spindle and kinetochore associated complex subunit 3; C13orf3, chromosome 13 open

reading frame 3; spindle and kinetochore-associated protein 3; MGC4832; RAMA1; C13orf3;

Entrez Gene ID	221150
mRNA Refseq	NM_001166017
Protein Refseq	NP_001159489
UniProt ID	Q8IX90
Chromosome Location	13q11
Function	protein binding;