



Human TCTN3 blocking peptide (CDBP2936)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-TCTN3 antibody
Antigen Description	This gene encodes a member of the tectonic gene family which functions in Hedgehog signal transduction and development of the neural tube. Mutations in this gene have been associated with Orofaciodigital Syndrome IV and Joubert Syndrom 18. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Sep 2012]
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	TCTN3 tectonic family member 3 [Homo sapiens]
Official Symbol	TCTN3

Synonyms	TCTN3; tectonic family member 3; C10orf61, chromosome 10 open reading frame 61; tectonic-3; DKFZP564D116; TECT3; tectonic 3; C10orf61; DKFZp564D116;
Entrez Gene ID	26123
mRNA Refseq	NM_001143973
Protein Refseq	NP_001137445
UniProt ID	Q6NUS6
Chromosome Location	10q23.33