



Human TMEM184A blocking peptide (CDBP2999)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-TMEM184A antibody
Antigen Description	TMEM184A (transmembrane protein 184A) is a protein-coding gene. An important paralog of this gene is TMEM184C.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	TMEM184A transmembrane protein 184A [Homo sapiens]
Official Symbol	TMEM184A
Synonyms	TMEM184A; transmembrane protein 184A; MGC9712; FLJ24011;

Entrez Gene ID	202915
mRNA Refseq	NM_001097620
Protein Refseq	NP_001091089
UniProt ID	Q6ZMB5
Chromosome Location	7p22.3