



Human TMEM214 blocking peptide (CDBP3003)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-TMEM214 antibody
Antigen Description	TMEM214 (transmembrane protein 214) is a protein-coding gene.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	TMEM214 transmembrane protein 214 [Homo sapiens]
Official Symbol	TMEM214
Synonyms	TMEM214; transmembrane protein 214; FLJ20254; FLJ39682;
Entrez Gene ID	54867

mRNA Refseq	NM_001083590
Protein Refseq	NP_001077059
UniProt ID	Q6NUQ4
Chromosome Location	2p23.3