



Human TMEM38A blocking peptide (CDBP3004)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-TMEM38A antibody
Antigen Description	TMEM38A (transmembrane protein 38A) is a protein-coding gene. GO annotations related to this gene include potassium channel activity. An important paralog of this gene is TMEM38B.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	TMEM38A transmembrane protein 38A [Homo sapiens]
Official Symbol	TMEM38A
Synonyms	TMEM38A; transmembrane protein 38A; trimeric intracellular cation channel type A; MGC3169;

TRIC A; TRICA; TRIC-A;

Entrez Gene ID	79041
mRNA Refseq	NM_024074
Protein Refseq	NP_076979
UniProt ID	Q9H6F2
Chromosome Location	19p13.11
Function	potassium channel activity;