



Human PAPD5 blocking peptide (CDBP3049)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-TRF4-2/PAPD5 antibody
Antigen Description	PAPD5 (PAP associated domain containing 5) is a protein-coding gene. GO annotations related to this gene include DNA-directed DNA polymerase activity and molecular_function. An important paralog of this gene is PAPD7.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	PAPD5 PAP associated domain containing 5 [Homo sapiens]
Official Symbol	PAPD5
Synonyms	PAPD5; PAP associated domain containing 5; PAP-associated domain-containing protein 5; TRF4 2; TUTase3; TUTase 3; terminal uridylyltransferase 3; topoisomerase-related function protein 4-2; TRF4-2; FLJ40270;
Entrez Gene ID	64282

mRNA Refseq	NM_001040284
Protein Refseq	NP_001035374
UniProt ID	Q8NDF8
Chromosome Location	16q12.1
Function	DNA binding; DNA-directed DNA polymerase activity; metal ion binding; molecular_function; nucleotidyltransferase activity; transferase activity;