



Human ZBED2 blocking peptide (CDBP3224)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-ZBED2 antibody
Antigen Description	ZBED2 (zinc finger, BED-type containing 2) is a protein-coding gene. Diseases associated with ZBED2 include gastric cancer. GO annotations related to this gene include DNA binding and metal ion binding. An important paralog of this gene is ZBED3.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	ZBED2 zinc finger, BED-type containing 2 [Homo sapiens]
Official Symbol	ZBED2
Synonyms	ZBED2; zinc finger, BED-type containing 2; zinc finger BED domain-containing protein 2; MGC10796; zinc finger, BED domain containing 2;

Entrez Gene ID	79413
mRNA Refseq	NM_024508
Protein Refseq	NP_078784
UniProt ID	Q9BTP6
Chromosome Location	3p13-q13.2
Function	DNA binding; metal ion binding;