



RapidGet™ Rat Immunoglobulin Isotyping Kit (CDXPH004)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Format	Cassette
Size	5T, 10T, 20T
Sample	hybridoma culture supernatant and purified monoclonal antibody
Intended Use	5 minute determination of 6 Rat Ig heavy chain isotypes and 2 light chain subclasses: IgG1, IgG2a, IgG2b, IgG2c, IgA, IgM, kappa, and lambda.
Assay Procedure	<ol style="list-style-type: none"> 1. Place the cassettes in horizontal plane. 2. Add 80-100 µL antibody solution directly to the Sample Well. Wait for 5 minutes. 3. A red test line at the location corresponding to antibody-specific isotype and a red control line at "C" position will appear within 1-5 minutes. <p>##NOTE##:</p> <ol style="list-style-type: none"> (1) Do not read after 20 minutes. (2) If the color of tissue culture medium causes high background, 40-60 µL PBS can be chased to clear the background. (3) Promptly return unused strips in their original bag and keep dry.