



Anti-C6ORF25 (full length) polyclonal antibody (CPBT-28818MH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Polyclonal antibody to Human C6ORF25.
Antigen Description	This gene is a member of the immunoglobulin (Ig) superfamily and is located in the major histocompatibility complex (MHC) class III region. The protein encoded by this gene is a glycosylated, plasma membrane-bound cell surface receptor, but soluble isoforms encoded by some transcript variants have been found in the endoplasmic reticulum and Golgi before being secreted. Multiple transcript variants encoding different isoforms have been found for this gene.
Immunogen	Full length human C6orf25 protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Isoforms A and B: Cell membrane; Single-pass type I membrane protein. Isoforms D and E: Endoplasmic reticulum. Golgi apparatus.
Format	Liquid
Size	50 µl
Buffer	Preservative: None Constituents: Whole serum

Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	C6orf25 chromosome 6 open reading frame 25 [Homo sapiens]
Official Symbol	C6ORF25
Synonyms	C6ORF25; chromosome 6 open reading frame 25; protein G6b; G6b; NG31; Chromosome 6 open reading frame 25; G6b; G6B protein; Immunoglobulin receptor; MGC142279; MGC142281; NG31; Protein G6b; immunoglobulin receptor;
Entrez Gene ID	80739
Protein Refseq	NP_079536
UniProt ID	B0V023
Chromosome Location	6p21
Function	heparin binding; receptor activity;