



Anti-TPBPA (aa 1-100) polyclonal antibody (CPBT-48527RM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Mouse TPBPA.
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 1 - 100 of Mouse Trophoblast specific protein alpha.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, ICC/IF, WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40Preservative:0.02% Sodium azideConstituent:PBSNote: Batches of this product that have a concentration
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	Tpbpa trophoblast specific protein alpha [Mus musculus]
Official Symbol	TPBPA
Synonyms	TPBPA; trophoblast specific protein alpha; trophoblast-specific protein alpha; Tb1; Tpbp; Tpbpa; Tpbpa trophoblast specific protein alpha; Trophoblast specific protein alpha; Trophoblast specific protein; clone 4311; Tb1; 4311; Tb-1; Tpbp; AV033951;
Entrez Gene ID	21984
Protein Refseq	NP_033437
UniProt ID	Q9CPR0
Function	molecular_function;