



Goat Anti-Human IgA polyclonal antibody [FITC] (CPBT-65467GA)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	CPBT-65467GA has cross absorption against human IgM and IgG. CPBT-65467GA might also cross react with IgA from other species.
Immunogen	Human IgA
Isotype	F(ab') ₂ IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Affinity purification
Conjugate	FITC
Applications	Suitable for use in IHC-F, FC, ELISA, and Stimulation. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	500 µg
Buffer	PBS
Preservative	0.1% Sodium Azide
Storage	2-8°C; Avoid exposure to light. This product is photosensitive and should be protected from

light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Ship

Wet ice

BACKGROUND

Keywords

IGHA2; immunoglobulin heavy constant alpha 2 (A2m marker); IGA;
