



Anti-Rat IgG polyclonal antibody [FITC] (CPBT-66021RA)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat IgG Flow Cytometry Use 50ul of the suggested working dilution to label 106 cells in 100ul.
Specificity	IgG
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	FITC
Applications	FC
Procedure	Secondary Antibodies
Size	1 mg
Preservative	0.09% Sodium Azide
Storage	in frost-free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Warnings	For research purposes only

GENE INFORMATION

Gene Name [Ighg Immunoglobulin heavy chain \(gamma polypeptide\) \[Rattus norvegicus \(Norway rat\) \]](#)

Official Symbol	Ighg
Synonyms	IGHG; Immunoglobulin heavy chain (gamma polypeptide); RGD1359539; IGG;
Entrez Gene ID	299354
UniProt ID	P20760
Chromosome Location	6q32
Pathway	Classical antibody-mediated complement activation; Complement cascade; Creation of C4 and C2 activators; FCGR activation; Fcgamma receptor (FCGR) dependent phagocytosis; Immune System; Initial triggering of complement; Innate Immune System; Regulation of actin dynamics for phagocytic cup formation; Role of phospholipids in phagocytosis;
Function	antigen binding