



# Anti-Rat IgG1 polyclonal antibody [AP] (CPBT-66039GR)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Goat anti Rat IgG1 has been shown to react with rat IgG1 by immunoelectrophoresis and ELISA. Goat anti Rat IgG1 may cross react with IgG1 from other species.
<b>Specificity</b>	IgG1
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	Rat
<b>Conjugate</b>	AP
<b>Applications</b>	ELISA; WB
<b>Format</b>	Purified IgG conjugated to Alkaline Phosphatase - liquid
<b>Size</b>	500 µg
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at +4°C. DO NOT FREEZE This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">Ighg Immunoglobulin heavy chain (gamma polypeptide) [ Rattus norvegicus (Norway rat) ]</a>
<b>Official Symbol</b>	Ighg

<b>Synonyms</b>	IGHG; Immunoglobulin heavy chain (gamma polypeptide); RGD1359539; IGG1;
<b>Entrez Gene ID</b>	<a href="#">299354</a>
<b>UniProt ID</b>	P20759
<b>Chromosome Location</b>	6q32
<b>Pathway</b>	Classical antibody-mediated complement activation; Complement cascade; Creation of C4 and C2 activators; FCGR activation; Fcgamma receptor (FCGR) dependent phagocytosis; Immune System; Initial triggering of complement;
<b>Function</b>	antigen binding