



Anti-Rabbit IgG Fc polyclonal antibody [FITC] (CPBT-67169GR)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat anti rabbit IgG antibody recognizes the Fc region of rabbit immunoglobulin G, confirmed by ELISA and immunoelectrophoresis. The antibody has not been species cross adsorbed and may cross react with IgG from other species. Goat anti rabbit IgG antibody has been used successfully as a secondary detection reagent for rabbit primary antibodies using immunohistochemical techniques.
-------------------------	--

Specificity	IgG
Isotype	IgG
Source/Host	Goat
Species Reactivity	Rabbit
Conjugate	FITC
Applications	FC
Format	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Size	1 mg
Preservative	0.09% Sodium Azide
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	LOC100009097 Ig gamma chain C region [Oryctolagus cuniculus (rabbit)]
Official Symbol	IGGC
Synonyms	LOC100009097; IGGC; Ig gamma chain C region; Ig gamma H-chain; immunoglobulin gamma1 constant region; IGG;
Entrez Gene ID	100009097
UniProt ID	P01870
Pathway	African trypanosomiasis;Allograft rejection; Amoebiasis; Asthma; Autoimmune thyroid disease;