



# Anti-Rabbit IgG polyclonal antibody [FITC] (CPBT-67174GR)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	This product recognizes rabbit Immunoglobulin G heavy and light chains. Purity Assay by IEP resulted in a single precipitin arc against anti-fluorescein, anti-goat serum, rabbit IgG and rabbit serum.
<b>Specificity</b>	IgG
<b>Immunogen</b>	Rabbit IgG whole molecule
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	Rabbit
<b>Conjugate</b>	FITC
<b>Applications</b>	FC; IF
<b>Size</b>	2 mg
<b>Preservative</b>	0.01% Sodium Azide
<b>Storage</b>	in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">LOC100009097 Ig gamma chain C region [ Oryctolagus cuniculus (rabbit) ]</a>
<b>Official Symbol</b>	IGGC
<b>Synonyms</b>	LOC100009097; IGGC; Ig gamma chain C region; Ig gamma H-chain; immunoglobulin gamma1 constant region; IGG;
<b>Entrez Gene ID</b>	<a href="#">100009097</a>
<b>UniProt ID</b>	P01839
<b>Pathway</b>	African trypanosomiasis; Allograft rejection; Amoebiasis; Asthma;