



Anti-Rabbit IgG polyclonal antibody [FITC] (CPBT-67188SR)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep anti Rabbit IgG recognizes rabbit IgG. Flow Cytometry Use 50ul of the suggested working dilution to label 106 cells in 100ul.
Specificity	IGGC
Immunogen	Rabbit IgG.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Rabbit
Conjugate	FITC
Applications	IHC-Fr; FC
Format	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Size	1 mg
Preservative	0.09% Sodium Azide
Storage	in frost-free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name [LOC100009097 Ig gamma chain C region \[Oryctolagus cuniculus \(rabbit\) \]](#)

Official Symbol	IGGC
Synonyms	LOC100009097; IGGC; Ig gamma chain C region; Ig gamma H-chain; immunoglobulin gamma1 constant region;
Entrez Gene ID	100009097
UniProt ID	P01870
Pathway	African trypanosomiasis;Allograft rejection; Amoebiasis; Asthma; Autoimmune thyroid disease;