



Anti-Rabbit IgG polyclonal antibody [R-PE] (CPBT-67189SR)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This product recognises rabbit IgG. Flow Cytometry Use 50ul of the suggested working dilution to label 106 cells in 100ul.
Specificity	IGGC
Immunogen	Rabbit IgG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Rabbit
Conjugate	PE
Applications	FC
Format	Purified IgG conjugated to R. Phycoerythrin (RPE) - Iyophilised
Size	100 tests
Preservative	0.09% Sodium Azide
Storage	Prior to reconstitution store at +4°C. Following reconstitution store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

45-1 Ramsey Road, Shirley, NY 11967, USA

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221

Gene Name	LOC100009097 Ig gamma chain C region [Oryctolagus cuniculus (rabbit)]
Official Symbol	IGGC
Synonyms	LOC100009097; IGGC; Ig gamma chain C region; Ig gamma H-chain; immunoglobulin gamma1 constant region;
Entrez Gene ID	100009097
UniProt ID	P01870
Pathway	African trypanosomiasis; Allograft rejection; Amoebiasis; Asthma; Autoimmune thyroid disease;