



Anti-Rabbit IgG polyclonal antibody [Dylight® 800] (CPBT-67193SR)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep anti Rabbit IgG antibody recognizes rabbit IgG. In ELISA techniques; 5% cross-reactivity is seen with human, mouse, rat or bovine serum. Flow Cytometry Use 50ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul.
Specificity	IGGC
Immunogen	Purified Rabbit IgG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Rabbit
Conjugate	Dylight 800
Applications	FC; IF; WB
Size	100 µg
Preservative	0.09% Sodium Azide
Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	LOC100009097 Ig gamma chain C region [Oryctolagus cuniculus (rabbit)]
Official Symbol	IGGC
Synonyms	LOC100009097; IGGC; Ig gamma chain C region; Ig gamma H-chain; immunoglobulin gamma1 constant region;
Entrez Gene ID	100009097
UniProt ID	P01870
Pathway	African trypanosomiasis;Allograft rejection; Amoebiasis; Asthma; Autoimmune thyroid disease;