



Anti-Rabbit IgG polyclonal antibody [HRP] (CPBT-67195SR)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This product recognises rabbit IgG as demonstrated by ELISA and double diffusion. Cross-reactivity against human, mouse, rat and bovine immunoglobulins is reduced by solid phase adsorption.
Specificity	IGGC
Immunogen	Purified rabbit IgG from normal rabbit serum
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Rabbit
Conjugate	HRP
Applications	IHC-Fr; ELISA; WB
Format	Purified IgG fraction conjugated to Horseradish Peroxidase (HRP) - liquid
Size	1 mg
Preservative	None
Storage	in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	LOC100009097 Ig gamma chain C region [Oryctolagus cuniculus (rabbit)]
Official Symbol	IGGC
Synonyms	LOC100009097; IGGC; Ig gamma chain C region; Ig gamma H-chain; immunoglobulin gamma1 constant region;
Entrez Gene ID	100009097
UniProt ID	P01870
Pathway	African trypanosomiasis;Allograft rejection; Amoebiasis; Asthma; Autoimmune thyroid disease;