



Anti-IgA polyclonal antibody [FITC] (CPBT-67965GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat anti Mouse IgA antibody recognizes Mouse IgA and has been cross absorbed against mouse IgM, IgG1, IgG2a, IgG2b and IgG3. Flow Cytometry Use 10ul of the suggested working dilution to label 106 cells in 100ul.
Specificity	IgA
Immunogen	Mouse IgA paraproteins
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse, Rat
Conjugate	FITC
Applications	IHC-Fr; FC; IF
Format	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Size	500 µg
Preservative	0.1% Sodium Azide
Storage	in frost free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	Igh-VJ558 immunoglobulin heavy chain (J558 family) [Mus musculus (house mouse)]
Official Symbol	Igh-VJ558
Synonyms	IGH-VJ558; immunoglobulin heavy chain (J558 family); AI893585; Igh-A (1g2); Vh186.2/Jh2; immunoglobulin heavy chain variable regionIGA;
Entrez Gene ID	16061
UniProt ID	P01878
Chromosome Location	12 F2; 12