



Anti-Mouse IgE polyclonal antibody [HRP] (CPBT-68020GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat anti mouse IgE polyclonal antibody recognizes mouse IgE and shows negligible cross - reactivity with other immunoglobulins.
Specificity	IgE
Immunogen	Native mouse IgE
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	HRP
Applications	ELISA
Format	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Size	500 µg
Preservative	None
Storage	in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name [Ighε Immunoglobulin heavy constant epsilon \[Mus musculus \(house mouse\) \]](#)

Official Symbol	Ighe
Synonyms	IGHE; Immunoglobulin heavy constant epsilon; Gm900; IGE;
Entrez Gene ID	380792
UniProt ID	P06336
Chromosome Location	12; 12 F1
Function	molecular_function