



## Anti-Mouse IgG1 polyclonal antibody [HRP] (CPBT-68024GM)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

| Product Overview   | Goat anti Mouse IgG1 antibody recognizes Mouse IgG1. This antibody has been cross absorbed against mouse IgM, IgG2a, IgG2b, IgG3 and IgA, pooled human sera and purified human paraproteins. Goat anti Mouse IgG1 antibody shows minimal cross-reactivity with human immunoglobulins. |
|--------------------|---|
| Specificity        | IGHG1   |
| Immunogen          | Mouse IgG1 paraproteins.  |
| Isotype            | IgG   |
| Source/Host        | Goat  |
| Species Reactivity | Mouse   |
| Conjugate          | HRP   |
| Applications       | IHC-Fr; ELISA; IHC-P; WB  |
| Format             | Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid  |
| Size               | 1 ml  |
| Buffer             | 50% Glycerol/50% Phosphate buffered saline, pH 7.4  |
| Preservative       | None  |
| Storage            | in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.  |

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 Fax: 1-631-938-8221

Email: info@creative-diagnostics.com

© Creative Diagnostics All Rights Reserved

## **GENE INFORMATION**

| Gene Name           | Ighg1 immunoglobulin heavy constant gamma 1 (G1m marker) [ Mus musculus (house mouse) |
|---------------------|---|
| Official Symbol     | IGHG1   |
| Synonyms            | IGHG1; immunoglobulin heavy constant gamma 1 (G1m marker); IgG1; Igh-4; VH7183;       |
| Entrez Gene ID      | <u>16017</u>  |
| UniProt ID          | P01868  |
| Chromosome Location | 12 F1; 12 62.1 cM   |
| Function            | antigen binding;  |