



Anti-Mouse IgG2a polyclonal antibody [FITC] (CPBT-68041GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	IGHG2A
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	FITC
Applications	IHC-Fr; FC
Size	400 µg
Preservative	None

GENE INFORMATION

Gene Name	Ighg2a immunoglobulin heavy constant gamma 2A [Mus musculus (house mouse)]
Official Symbol	Ighg2a
Synonyms	Ighg2a; immunoglobulin heavy constant gamma 2A; Igh-1; Igh-1a; 1810060O09Rik; Ig heavy chain Mem5; MLT2-23 heavy chain; immunoglobulin gamma-2a heavy chain; immunoglobulin heavy chain 1 (serum IgG2a); immunoglobulin heavy chain 1a (serum IgG2a)
Entrez Gene ID	380793
Chromosome Location	12 F2; 12

Functionantigen binding
