



Anti-Mouse IgM polyclonal antibody [HRP] (CPBT-68045GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat anti Mouse IgM antibody recognizes murine IgM and is cross-adsorbed to reduce reactivity with other mouse heavy or light chains.
Specificity	IGHM
Immunogen	Mouse IgM mu heavy chain.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	HRP
Applications	IHC-Fr; ELISA; IHC-P
Format	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Size	1 mg
Preservative	None
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name [Ighm immunoglobulin heavy constant mu \[Mus musculus \(house mouse\) \]](#)

Official Symbol	Ighm
Synonyms	IGHM; immunoglobulin heavy constant mu; Igm; muH; Igh6; Igh-6; Igh-M; AI326478; Ig mu; immunoglobulin heavy chain 6 (heavy chain of IgM); mu-immunoglobulin;
Entrez Gene ID	16019
UniProt ID	P01873
Chromosome Location	12 62.1 cM; 12 F1-2
Pathway	Adaptive Immune System; Antigen activates B Cell Receptor (BCR) leading to generation of second messengers; Immune System; Signaling by the B Cell Receptor (BCR);
Function	antigen binding; protein binding; transmembrane signaling receptor activity