



Anti-Mouse IgG Fc polyclonal antibody [FITC] (CPBT-68072GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	IgG
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	FITC
Applications	IHC-Fr; FC
Size	1 mg
Preservative	None

GENE INFORMATION

Gene Name	IGKC immunoglobulin kappa constant [Homo sapiens (human)]
Official Symbol	IGKC
Synonyms	IGKC; immunoglobulin kappa constant; Km; HCAK1; IGKCD; IGG;
Entrez Gene ID	<u>3514</u>
UniProt ID	P01834
Chromosome Location	2p12

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221

Pathway

Adaptive Immune System; Antigen activates B Cell Receptor (BCR) leading to generation of second messengers; Binding and Uptake of Ligands by Scavenger Receptors; Classical antibody-mediated complement activation; Complement cascade; Creation of C4 and C2 activators; FCERI mediated Ca+2 mobilization; FCERI mediated MAPK activation;

Function

antigen binding;

Tel: 1-631-624-4882 Fax: 1-631-938-8221

2/2