



Anti-IgG polyclonal antibody [HRP] (CPBT-68077RA)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse IgG (all subclasses). Some cross reactivity with IgM will be expected. Cross reactivity is expected with rat IgG. Cross reactivity with human serum has been minimised by solid phase adsorption.
Specificity	IGHG1
Immunogen	Purified whole mouse immunoglobulin G
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat
Conjugate	HRP
Applications	IHC-Fr; ELISA; IHC-P; IHC-R; WB
Format	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Size	1 mg
Preservative	See individual product datasheet
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	Ighg1 immunoglobulin heavy constant gamma 1 (G1m marker) [Mus musculus (house mouse)]
Official Symbol	IGHG1
Synonyms	IGHG1; immunoglobulin heavy constant gamma 1 (G1m marker); IgG1; Igh-4; VH7183;
Entrez Gene ID	16017
UniProt ID	P01863
Chromosome Location	12 F1; 12 62.1 cM
Function	antigen binding;