



Anti-Mouse IgG polyclonal antibody [HRP] (CPBT-68087GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This product recognises mouse IgG. Cross-reactivity with rat immunoglobulins has been reduced by solid phase adsorption.
Specificity	IGHG1
Immunogen	Mouse IgG
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	HRP
Applications	IHC-Fr; ELISA; IHC-P
Format	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Size	500 µg
Preservative	See individual product datasheet
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	Ighg1 immunoglobulin heavy constant gamma 1 (G1m marker) [Mus musculus (house mouse)]
------------------	---

Official Symbol	IGHG1
Synonyms	IGHG1; immunoglobulin heavy constant gamma 1 (G1m marker); IgG1; Igh-4; VH7183;
Entrez Gene ID	16017
UniProt ID	P01863
Chromosome Location	12 F1; 12 62.1 cM
Function	antigen binding;